



InTray™ Colorex™ MRSA

Instructions for Use

20-3301 InTray™ Colorex™ MRSA, 2"tray, 20 trays/box

20-3307 InTray™ Colorex™ MRSA, 2"tray, 5 trays/box

INTENDED USE

Colorex™ MRSA is a selective and differential chromogenic medium, containing selective agents that exhibit high sensitivity and specificity for the isolation of Methicillin Resistant *Staphylococcus aureus*. Intended for use with a direct streaking technique to differentiate Methicillin Resistant *Staphylococcus aureus*, the test can be performed with common swab samples composed of mixed populations of bacteria, e.g., nasal, perineal, throat, rectal specimens, etc. Colorex™ MRSA is intended for use in the identification of colonization with Methicillin Resistant *Staphylococcus aureus* to aid in the prevention and control of MRSA in healthcare settings. Colorex™ MRSA is not intended to diagnose infections by Methicillin Resistant *Staphylococcus aureus*, guide or monitor treatment for infections, or provide susceptibility results for Methicillin Resistant *Staphylococcus aureus*. Further confirmatory identification as Methicillin Resistant *Staphylococcus aureus*, e.g., by biochemical tests such as latex agglutination, *mecA* gene PCR, and/or disk diffusion method for oxacillin and cefoxitin, is required (1).

DESCRIPTION AND PRINCIPLE OF USE

Methicillin Resistant *Staphylococcus aureus* exhibits resistance to a large panel of antibiotics, including beta-lactam antibiotics. MRSA is one of the leading causes of nosocomial infections. Sources are either endogenous (the patient) or through cross contamination (environmental or person to person contact). Pre-admission and wide-scale screening for MRSA has proved to be an effective method for reducing hospital burden of MRSA-colonized patients while reducing response time and laboratory workload.

STORAGE

Upon receipt, store InTray™ Colorex™ MRSA under refrigeration (2-8°C). Medium can be kept for one day at ambient temperature. Protect media from exposure to light, excessive heat, moisture and freezing. Do not open until ready to use. Do not use if the medium shows signs of deterioration, shrinking, cracking, discoloration or contamination.

INOCULATION PROCEDURE

Allow the InTray™ to warm to room temperature. Lift the lower right corner of the flexible InTray™ label until the protective seal is completely visible. Remove the paper-foil seal by pulling the tab (Fig. 1). **Discard** the seal. **DO NOT REMOVE OR ALTER THE WHITE FILTER STRIP OVER THE VENT HOLE.**




Figure 1




Figure 2



Figure 3

Streak sample onto the agar surface for isolation (Fig. 2). Reseal the InTray™ label to the plastic tray body. **Press all around the perimeter of the InTray™ to ensure a complete seal** (Fig. 3). Immediately label the InTray™ with sample information and date. **DO NOT COVER THE VIEWING WINDOW.**

CULTURE AND RESULTS

Incubate at 37°±2°C for 18-24 hours under ambient atmosphere. Colonies of Methicillin Resistant *Staphylococcus aureus* appear rose to mauve. Methicillin Susceptible *Staphylococcus aureus* is inhibited. Other non-resistant bacteria or yeast are blue, colorless or inhibited.



LIMITATIONS/PRECAUTIONS

For *in vitro* diagnostic use only by compliant and trained professional laboratory personell. InTray™ Colorex™ MRSA is not intended to diagnose infections by Methicillin Resistant *Staphylococcus aureus*, guide or monitor treatment for infections, or provide susceptibility results for Methicillin Resistant *Staphylococcus aureus*. Further confirmatory identification as Methicillin Resistant *Staphylococcus aureus*, e.g., by biochemical tests such as latex agglutination, mecA gene PCR, and/or disk diffusion method for oxacillin and ceftioxin, is required (1). Do not use any InTrays™ that are beyond its expiry date. Because of the potential for containing infectious materials, the InTray™ must be destroyed by autoclaving at 121°C for 20 minutes.

INTERPRETATION

Organism	Colony Appearance
MRSA strain <i>S. aureus</i>	Rose to Pink-Mauve with matte halo
MSSA strain <i>S. aureus</i>	Inhibited
Other gram (-), (+) bacteria and yeast	Inhibited/blue/colorless



REAGENTS

Colorex™ MRSA contains agar, peptone nutrients, salts, antimicrobial selective compounds and chromogenic additives.

QUALITY CONTROL

All Colorex™ MRSA Agar products are performance verified with the following ATCC® microbe strains. Product performance is also verified periodically throughout the marked shelf life of each lot.

Organism	ATCC®	Colony Appearance
MRSA Strain: <i>S. aureus</i>	BAA-1720	Rose to Pink-Mauve with matte halo
MSSA strains: <i>S. aureus</i>	25923	Inhibited
Other gram (-), (+): <i>P. aeruginosa</i>	9027	“
<i>E. faecalis</i>	29212	“
<i>E. coli</i>	25922	“
Yeast: <i>C. albicans</i>	60193	“

SYMBOL KEY			
Symbol	Used For	Symbol	Used For
	Batch code		Temperature limitation
	Date of manufacture		Catalog number
	Use by YYY-MM-DD or YYYY-MM		Caution, consult accompanying documents
	Manufacturer		Authorized representative in the European Community
	In vitro diagnostic medical device		in European community

1. Goodwin KD and M Pobuda (2009) Performance of CHROMagar™ Staph aureus and CHROMagar™ MRSA for detection of *Staphylococcus aureus* in seawater and beach sand – Comparison of culture, agglutination, and molecular analyses. *Water Research* 43 (4802–4811).

Biomed Diagnostics, Inc.
 PO Box 2366 • 1388 Antelope Road
 White City, OR 97503
 (800) 964-6466 • (541) 830-3000
 www.biomeddiagnostics.com
 info@biomeddiagnostics.com

COLOREX™ is a trademark of Dr. A. Rambach, France
 ATCC® is a registered trademark of the American Type Culture Collection