

# InTray™ BG Agar

## InTray™ BG Agar + Novobiocin

Brilliant Green Agar and  
Brilliant Green Agar with Novobiocin antibiotic

40-1007 (5 pack); 40-1001 (20 pack)

For *In Vitro* Diagnostic Use Only



### INTENDED USE:

Brilliant Green Agar is a highly selective medium for the isolation of *Salmonella* other than *S. typhi* from feces and other materials. Brilliant Green Agar with addition of Novobiocin confers an additional measure of selectivity by inhibiting contaminant organisms commonly seen on agar media when isolating *Salmonella*<sup>1</sup>.

### DESCRIPTION AND PRINCIPLE OF USE:

Brilliant Green Agar was first described by Kirstensen in 1925<sup>2</sup> and later improved by Kauffmann in 1935<sup>3</sup>. It is a highly selective medium for the recovery of salmonellae except for the typhoid and paratyphoid bacilli strains<sup>4</sup>. Brilliant Green dye inhibits gram-positive bacteria and a majority of gram-negative bacilli. Phenol red serves as a pH indicator, yielding a yellow color as a result of acid formation from lactose and sucrose fermentation.

### INOCULATION PROCEDURE:

Pull back the lower right corner adjacent to the clear window of the InTray™ label until the protective seal is completely visible. Remove the seal by pulling the tab (Fig. 1). **Discard** the seal. **DO NOT REMOVE OR ALTER THE WHITE FILTER STRIP OVER THE VENT HOLE!**




Figure 1

Inoculate the specimen on the surface of the medium. A sterile inoculating loop that has been moistened by touching the surface of the medium may be used for inoculation of solids or scrapings (Fig. 2).




Figure 2

Reseal the InTray™ by pressing together the edges of the label against the plastic tray. **Press all around the InTray™ to insure a complete seal** (Fig. 3). Immediately label the InTray™ with patient or sample information and date. **DO NOT COVER THE VIEWING WINDOW.**



Figure 3

### CULTURE AND RESULTS:

Use standard procedures to obtain isolated colonies from specimens. A less-selective or non-selective medium should also be streaked to increase the chance of recovery when the population of gram-negative bacteria in the sample may be low.

**Incubate aerobically at 35 ± 2°C for 18-24 hours; If negative after 24 hours, re-incubate an additional 24 hours.**

### CULTURE RESPONSE:

Organism	ATCC™	Recovery	Colony Color
<i>Escherichia coli</i>	25922	Poor	Yellow-Green
<i>Salmonella choleraesuis</i> (Enteridis serotype)	13076	Good	Red
<i>Salmonella choleraesuis</i> (Typhi serotype)	19430	None to Poor	Red
<i>Salmonella choleraesuis</i> (Typhimurium serotype)	14028	Good	Red
<i>Staphylococcus aureus</i>	25923	Marked inhibition	-

### STORAGE:

Upon receipt, store InTray™ BG and BG+Novobiocin under refrigeration (2-8°C). Avoid freezing or prolonged storage at temperatures above 40°C. Do not use InTray™ BG or BG+Novobiocin if the medium shows signs of deterioration or contamination. Expiration is 12 months past the date of manufacture.

### PRECAUTIONS:

For *in vitro* diagnostic use only. Once the InTray™ has been inoculated and resealed, open only in a biological safety cabinet. Because of the potential for containing infectious materials, the InTray™ must be destroyed by autoclave or other equivalent means of disposal.

SYMBOL KEY			
Symbol	Used For	Symbol	Used For
	Batch code		Temperature limitation
	Date of manufacture		Catalog number
	Use by YYY-MM-DD or YYYY-MM		Caution, consult accompanying documents
	Manufacturer		Authorized representative in the European Community
	In vitro diagnostic medical device		in European community

### REFERENCES:

1. Miller and Tate. 1990. Maryland Poultryman..
2. Kirstensen *et al.* 1925. *Br. J. Exp. Pathol.*
3. Kauffmann. 1935. *Z. Hyg. Infektionskr.*
4. Sack *et al.* *Cumitech 12*, American Society for Microbiology, 1980.

BioMed Diagnostics, Inc.

PO Box 2366 • 1388 Antelope Road • White City, OR 97503

(800) 964-6466 / (541) 830-3000

www.biomeddiagnostics.com