INTRODUCTION

INTENDED USE

InTray® SAB-FungID contains Sabouraud's dextrose agar with chloramphenicol and cycloheximide, a selective medium used to aid in the detection of dermatophyte fungi from clinical specimens with mixed microbiota.

DESCRIPTION & PRINCIPLE OF USE

The InTray® SAB-FungID is a dynamic system with built-in components and features that are designed for user compatibility and ease of dermatophyte fungi detection (i.e., 10x Obj. microscopy direct from the tray). The medium allows for the growth & observation of distinct colony morphology and color (i.e., pigments) of dermatophyte fungi, while inhibiting most gram-positive bacteria, gram-negative bacteria, yeast and saprophytic fungi. Dermatophytes are fungi in the genera *Microsporum, Tricophyton* and *Epidermophyton*. They are capable of metabolizing keratin found in skin, hair and nails of living hosts. The fungi characteristically may invade the cutaneous tissue of the living host but rarely penetrate the subcutaneous tissue. Tinea and ringworm are two terms commonly used to describe dermatophytes.

REAGENTS & APPEARANCE

Sabouraud's dextrose agar appears transparent with a light amber hue & contains peptic/casein digest, dextrose, chloramphenicol (0.050g/L) & cycloheximide (0.40g/L) with a final pH of 5.6 ± 0.2 @ 25°C.

STORAGE & SHELF LIFE

On receipt, store trays under refrigeration or at room temperature (2- 25° C) in the dark. Avoid freezing or prolonged storage at temperatures greater than 40°C. Do not use trays if the medium shows signs of deterioration or contamination.

KEY NOTES REGARDING SPECIMEN COLLECTION

Specimen collection poses a major uncertainty in using this device.

NAILS - Often collecting viable material from infected nails is difficult because the living organisms are well under the nail itself. For best results, cut nails into small pieces.

HAIR - samples should be grasped at the uninfected end and several (3-6) small pieces, about 2 cm long, should be cut from the infected portion for inoculation onto the medium.

SKIN - scrapings should be taken with an inoculation tool that has been moistened with the medium or a sharp blade from the outer ridge of an active lesion. Vesicular fluid is unacceptable for dermatophyte culture. If the infected area is vesiculated, skin scrapings should be taken from the surface.

INSTRUCTIONS FOR USE

SAMPLE COLLECTION

Gather materials needed for the test:

- InTray® SAB-FungID test(s)
- Disposable gloves and sterile inoculation tool (e.g., cotton swab/forceps/scalpel blade)
- Laboratory incubator

<u>Prepare Sample</u>: Use aseptic technique during specimen collection and handling. Remove any soap residue from the sampling area. Clean the area with 70% alcohol and permit to air dry.

Collect Sample: InTray® SAB-FungID is designed for culturing hair, skin and nail samples (i.e., cuttings/scrapings). All specimens should be handled according to CDC infectious materials isolation guidelines;

https://www.cdc.gov/infectioncontrol/guidelines/isolation/index.html

INOCULATION

Allow tray(s) to warm to room temperature before use. Immediately label the tray with patient/sample information and date. Pull back the lower right corner adjacent to the clear window of the tray label until the protective seal is completely visible. Remove the seal by pulling the tab. **Discard** the seal. **DO NOT REMOVE OR ALTER THE WHITE FILTER STRIP OVER THE VENT HOLE!**



Inoculate the specimen on the **center** surface of the medium. A sterile inoculating loop that has been moistened by touching the surface of the medium may be used for inoculation of solids or scrapings.



Re-seal all around the tray to ensure a complete seal by pressing the edges of the label against the plastic tray. <u>DO NOT COVER THE</u> <u>VIEWING WINDOW.</u> Complete re-seal prevents dehydration!



INCUBATION

Incubate inoculated trays in a dark humidified environment for up to 21 days at 25°C-30°C. Observe the trays daily through the clear viewing window.

PRECAUTIONS, SAFETY & DISPOSAL

IVD For professional laboratory use only. Inoculated trays have the potential of containing live infectious materials, and they must be handled and destroyed in accordance with BSL-2 requirements. Destroy trays by autoclaving at 121°C for 20 minutes or through other suitable means of sterilization.

READING THE RESULTS

MACRO & MICROSCOPIC EVALUATION

Observe the medium for fungal colony growth, texture, morphology and pigmentation characteristic of dermatophytes and saprophytes. Without opening the InTray® SAB-FungID, place the tray under a 10x objective microscope lens to view distinct fungal structures (i.e., hyphae, micro/macro-conidia). See the Quick Identification Chart on the following page for reference. Trays to be used with 10x Obj. ONLY!

Positives: If, within 1 - 21 days, the medium shows characteristic dermatophyte fungi colony growth, texture, morphology, pigmentation and microscopic structure, the InTray® SAB-FungID is presumptive positive.

Negatives: Trays that show no colony growth after 21 days from the date of inoculation are presumptive negative.

LIMITATIONS

If a fungal infection is strongly suspected and the culture test result is negative, it may be appropriate to retest giving more care to specimen collection.

- Optimally, simultaneously inoculate samples on both InTray® DM-FungID & InTray® SAB without antimicrobials (e.g., for some pathogenic fungi inhibited by antimicrobials)
- Transfer of fungal growth to slide examination: Stained slide preparations can be made as appropriate based on 10x Obj. examination of InTray® SAB-FungID. InTray® Potato Dextrose Agar (PDA) can be used to induce sporulation.
- Final culture identification is to be made by laboratory professionals only, and may require further biochemical/culture testing.

QUALITY CONTROL

This product has been tested and meets the CLSI Approved Standard for commercially prepared media (M22-A3). At the time of manufacture, quality control testing is performed on each lot of InTray® SAB-FungID. The ability of the media to support growth and morphology is summarized on the Biomed Certificate of Analysis (CoA).

Table 1: Recommended Strains for QC Testing SAB-FungID

Test Strain	ATCC ^{тм} Number	Expected Results
T. mentagrophytes	9533	Growth
T. rubrum	28188	Growth
C. albicans	60193	Growth
A. brasiliensis	16404	Significant Inhibition
S. aureus	25923	Significant Inhibition
E. coli	25922	Significant Inhibition

OUICK IDENTIFICATION CHART (DERMATOPHYTES)

Epidermophyton floccosum

Colony appearance; green, yellow edges Reverse appearance; yellowish brown Septate hyphae. Macroconidia: (7-12 x 20-40µm) smooth, thick and thin-walled, club shaped with rounded ends, two



to six cells, singly or clusters. Microconidia: none.

Microsporum canis



Reverse appearance; yellow-orange Septate hyphae. Macroconidia: (10-25 x 35-110µm) numerous, long, spindle-shaped, rough, thick walls are apparent at the knob-like ends. Microconidia: few in number, smooth and club-shaped.

Colony appearance; white, orange-yellow

Microsporum gallinae



Reverse appearance; red Septate hyphae. Macroconidia are clavate to cigar shaped (6-8 x 15-50µm) and 2 to 10 celled. Microconidia are unicellular and ovoid to pyriform in shape.

Colony appearance; white, pink



Colony appearance; brown cinnamon Reverse appearance; light tan Septate hyphae. Macroconidia: (8-16 x 22-60µm) large numbers, symmetrical, rough and relatively thin-walled with rounded ends, not pointed like M. canis.

Microconidia: usually resent, club-shaped.

Microsporum nanum

Colony appearance; white/yellow/brown Reverse appearance; orange to red-brown

Reverse appearance; yellow to red



Septate hyphae. Macroconidia (4-8 x 12-18µm) are rough,



fairly thin-walled (as in M. gypseum), egg-shaped with a truncated base, having usually two cells. Microconidia are club shaped and smooth-walled. Abundance may vary. Colony appearance; white, tan

Trichophyton equinum

Abundant microconidia which may be clavate to pyriform 0 and sessile or spherical and stalked are formed laterally 00 along the hyphae. Macroconidia are rarely produced, but 200 variable size.

when present are clavate, smooth, thin walled and of Colony appearance; white to yellow

Trichophyton mentagrophytes



Reverse appearance; reddish brown Septate hyphae. Macroconidia: (4-8 x 20-50µm) occasionally present, cigar-shaped, thin-walled, narrow attachments to septate hyphae, 1 - 6 cells, found in young cultures 5 - 10 days old. Microconidia: usually present in powderv cultures, verv round, clustered on branched

Colony appearance; white, pink

Colony appearance; white, tan, yellow

Reverse appearance; yellow, tan

conidiophores; in fluffy cultures, smaller, fewer, teardrop shaped and easily confused with those of T. rubrum.

Trichophyton rubrum



Reverse appearance; yellow to wine red Septate hyphae. Macroconidia: (4-6 x 15-30µm) abundant, rare or absent, however can be long, narrow, thin-walled, parallel sides, 2 - 8 cells, may form on the ends singly or in groups. Microconidia: (2-3 x 3-5µm) lateral, teardrop shaped, form on macroconidia.

Trichophyton tonsurans

Septate hyphae with variable shaped microconidia all along the hyphae or on short conidiophores that are perpendicular to the parent hyphae. Microconidia are usually teardrop or club-shaped, but may elongate or enlarge to round "balloon" forms. Macroconidia are rare, irregular in form

QUICK IDENTIFICATION CHART (SAPROPHYTES/ CONTAMINANTS)

Alternaria sp. (Saprophyte)



Hyphae are septate and dark. Conidiophores are septate, variable in length and sometimes branched. Conidia are large (7-10 by 23-24µm), brown, have both transverse and longitudinal stations, and found singly or in chains. They are usually round at the end nearest the conidiophore, producing a club-like shape.

Aspergillus sp. (Saprophyte)



Microscopic morphology - Septate hypae (2.5 - 8µm in diameter); unbranched condidiophore arises from a specialized foot cell. The conidiophore is enlarged at the tip, forming a swollen vesicle completely or partially covered with flaskshaped phialides. The phialides produce chains of mostly round, sometimes rough conidia (2-5µm in diameter)

Candida albicans & Candida tropicalis



Microscopic examination shows clusters of bubbles in the media. Yeast is a common contaminant and there are inhibitors in the media to prevent most yeast species from growing. However, yeast species can be pathogenic depending upon the patient's overall condition and medication history.

Penicillium sp. (Saprophyte)



Microscopic morphology - Septate hyphae (1.5-5µm in diameter) with branched conidiophores that have secondary branches known at metulae. On the metulae, are flask-shaped phialides that bear unbranched chains of smooth or rough conidia (2.5-5µm in diameter). The entire structure forms the characteristic "penicillus" or "brush" appearance.

REFERENCES

1) Sabouraud R (1892). Ann Dermatol. Syphil. 3:1061.

2) Tille, et al. (2014) Bailey & Scott's Diagnostic Microbiology, Elsevier: p. 450. 3) Kwon-Chung, KJ and Bennett, JE (1992) MedicalMycology, Lea and Febiger: Philadelphia.

4) Murray, PR, Baron, ET, Pfaller, MA, Tenover, FC, Yolkem, RH, (1995) Manual of Clinical Micro- biology 6th ed., American Society for Microbiology: Washington, D.C., pp. 709-722.

5) Larone, DH (1995) Medically Important Fungi: A Guide to Identification, 2nd ed., American Society for Microbiology: Washington, D.C.

InTray® SAB-FungID

Sabouraud's Dextrose Agar with Cycloheximide & Chloramphenicol

Catalog No. 18-1207	5 Test Kit
Catalog No. 18-1201	20 Test Kit

A SELECTIVE CULTURE SYSTEM FOR THE IDENTIFICATION OF Dermatophyte fungi

For In Vitro Diagnostic Use





Manufactured by:

Biomed Diagnostics, Inc.

PO Box 2366 • White City, OR 97503 tel. (800)-964-6466 • fax. (541) 830-3001 info@biomeddiagnostics.com www.biomeddiagnostics.com

100-539 InTray® SAB Insert Rev. New (10/2017)