



Certificate of Analysis
InTray™ COLOREX™ Yeast Test

Product Number: 10-6101, 10-6107

Lot Number: 5CA036

Expiration Date: 09 MAR 2016

Test	Specification	Result
Appearance	Clear, colorless to light amber	Clear, colorless
Sterility	No growth on negative control	No Growth
pH	6.1 ± 0.2	6.02
Microbiology:		
<i>Candida glabrata</i>	Dark Rose	Rose
<i>Candida krusei</i>	Light pink / White edges	Light pink / white edges
<i>Candida albicans</i>	Green to Blue-Green	Green / blue
<i>Candida tropicalis</i>	Blue-Violet	Purple

QC Acceptance Date: 11 MAR 2015

Signature: 