



Certificate of Analysis
InPouch™ TV *Trichomonas vaginalis* Test

Product Number 10-2010 10-2003

Lot Number: 5DA045A

Expiration Date: 09 APR 2016

| Test | Specification | Result |
|---------------|--|---------|
| Appearance | Clear, Straw to Amber | Amber |
| pH Test | 6.05 to 6.15 | 6.10 |
| Sterility | Aseptic Medium Production and Dispensing | Sterile |
| Doubling Time | 3 – 8 Hours | 4.53 |

QC Acceptance Date: 20 APR 2015

Signature: 