



Certificate of Analysis
InPouch™ TV *Trichomonas vaginalis* Test

Product Number 10-2010 10-2003

Lot Number: **5EA045B**

Expiration Date: **09 APR 2016**

Test	Specification	Result
Appearance	Clear, Straw to Amber	Amber
pH Test	6.05 to 6.15	6.10
Sterility	Aseptic Medium Production and Dispensing	Sterile
Doubling Time	3 – 8 Hours	4.53

QC Acceptance Date: 20 APR 2015

Signature: 