



**Certificate of Analysis**  
**InPouch™ TV *Trichomonas vaginalis* Test**  
Product Number 10-2010 10-2003

Lot Number: 5FA059A

Expiration Date: 07 MAY 2016

Test	Specification	Result
Appearance	Clear, Straw to Amber	Clear, amber
pH Test	6.05 to 6.15	6.07
Sterility	Aseptic Medium Production and Dispensing	Sterile
Doubling Time	3 – 8 Hours	4.27

**QC Acceptance Date: 17 May 2015**

**Signature:** 