



Certificate of Analysis

InTray™ GC Agar

Product Number: 10-8001, 10-8007

Lot Number: **5FA091**

Expiration Date: **23 June 2016**

Test	Specification	Results
Appearance	Chocolate brown, opaque	Chocolate brown
Sterility	No growth on negative control	No growth
pH	7.2 ± 0.2	7.14
Microbiology (ATCC®):	Colony Appearance	
<i>N. gonorrhoeae</i> 43069	Colony growth	Growth
<i>N. meningitidis</i> 13090	Colony growth	Growth
<i>N. sicca</i> 9913	No Growth / Significant inhibition	No growth
<i>E. coli</i> 25922	No Growth / Significant inhibition	No growth
<i>P. mirabilis</i> 43071	No Growth / Significant inhibition	Significant inhibition
<i>S. epidermidis</i> 12228	No Growth / Significant inhibition	No growth
<i>C. albicans</i> 60193	No Growth / Significant inhibition	No growth

QC Acceptance Date: 29 June 2015

Signature: 