



Certificate of Analysis
InTray™ COLOREX™ Yeast Test

Product Number: 10-6101, 10-6107

Lot Number: **5GA092**

Expiration Date: **06 Jul 2016**

Test	Specification	Result
Appearance	Clear, colorless to light amber	Transparent, light amber
Sterility	No growth on negative control	No growth
pH	6.1 ± 0.2	6.08
Microbiology:		
<i>Candida glabrata</i>	Dark Rose	Rose
<i>Candida krusei</i>	Light pink / White edges	Light pink/white edges
<i>Candida albicans</i>	Green to Blue-Green	Green
<i>Candida tropicalis</i>	Blue-Violet	Blue/purple

QC Acceptance Date: 08 Jul 2015

Signature: 