



Certificate of Analysis
InPouch™ TF *Tritrichomonas foetus* Test

Product Number 11-1001, 11-1003

Lot Number: 5HA102A

Expiration Date: 23 July 2016

Test	Specification	Result
Appearance	Clear, Straw to Amber	Clear, amber
pH Test	6.6 to 6.8	6.74
Sterility	Aseptic Medium Production and Dispensing	No Growth
Doubling Time	3 – 8 Hours	4.17

QC Acceptance Date: 03 August 2015

Signature: 