


Learn more at biomeddiagnostics.com

- Collect, transport, incubate and result in one device
- View directly on the microscope stage
- Rapid TAT
- No re-plating required
- FDA Registration and Euro CE Mark*
[*on select products]
- Extended 6 - 27 month shelf life
[varies depending on product]
- Reduces contamination risk and lab materials



Resealable Label

Prevents media dehydration
Clear anti-fog view window

Media Well

Large 2-inch surface area
Produce distinctive colony growth

Durable Stackable Tray

Minimal Footprint