



# In partnership with Rhino Diagnostics, Biomed now offers **iClean<sup>®</sup> Nasopharyngeal Swabs**

Individually-sealed, sterile, flocked, nasopharyngeal swab with plastic, breakable shafts.  
For collecting SARS-CoV-2 and other clinical specimens.

## Product Principles

- Individually-sealed, sterile, flocked, nasopharyngeal swab with plastic, breakable shaft
- Proprietary nylon tip flocking process optimizes specimens collection and instant, spontaneous elution in media with up to 95% release
- No-slip, frosted, flexible plastic handle with molded breakpoint, ergonomic & anatomic designs for optimal patient comfort and efficient sample collection
- Compatible with most diagnostic molecular and antigen tests

## FDA Status

Manufacturer and product registered, listed, and FDA compliant

## Manufacturer

Chenyang Global is a leading medical consumables manufacturer with CE, FDA, ISO9001, ISO13485, and GMP Class 100,000 Cleanroom certifications. Supplier of COVID-19 supplies to numerous U.S. and European government agencies, healthcare providers, and World Health Organization

## Intended Use

For the collection of nasopharyngeal specimens, including SARS-CoV-2.

## Tip Type

100% medical grade nylon

## Tip Length

20 mm

## Tip Diameter

1.8 ± 0.2 mm



Learn more at  
[biomeddiagnostics.com](https://www.biomeddiagnostics.com)