

InPouch[®] TF-Bovine

The Gold Standard Trich test, InPouch TF is a nutrient both media, allows you to inoculate, culture and detect in one device. With a highly selective proprietary medium for *Tritrichomonas foetus*, InPouch's specificity is increased by inhibiting the growth of yeasts, molds and other bacteria.

PRINCIPLES

- Direct inoculation to media at the POC
- Easy detection and identification

COST SAVINGS

- Conveniently collect, culture and result in one device
- One-year shelf life at room temperature (18 - 25°C) from date of manufacture
- Oxygen resistant plastic pouch seals for transport and incubation

SPECIFICITY

- Specially formulated medium is selective for the growth of *T. foetus*, while inhibiting fungi and bacteria

SPECIMEN

Vaginal swab | Distal penile swab

RESULTS

- Evaluate 2-4 days for confirmed presence of Trichomonads.
- If no organisms or motions are detected, incubate the sample at 35 - 37°C.
(8% of positive results will occur within 5 days.)

AVAILABILITY

10 pack | 12-011-002
100 pack | 12-011-004

Also available in a PCR compatible transit tube:

10 pack | 12-012-001
50 pack | 12-012-002



Learn more at
biomeddiagnostics.com

or call 866.964.6466
to speak to a representative