

InPouch[®] TF-Feline

A gold standard feline trich test, InPouch TF-Feline allows you to inoculate, culture and detect in one device.

Created for use in your own clinic, InPouch is a sealed oxygen resistant, 2-chamber plastic pouch. This two-compartment system allows direct (wet mount) observation of newly inoculated specimens in the upper chamber before expressing it into the lower chamber for culture at 15 - 37°C.

Quality tested for sensitivity, InPouch's media can potentially result in a presumptive positive with an inoculum containing as little as one organism — all housed within a damage resistant form that maintains specimen viability during transport or incubation.

PRINCIPALS	- Direct inoculation to media at the POC - Easy detection and identification
COST SAVINGS	- Conveniently collect, culture and result - PCR compatible - One-year shelf life at room temperature (18 - 25°C) from date of manufacture - Oxygen resistant plastic pouch seals for transport and incubation
SPECIFICITY	- Specially formulated medium is selective for the growth of <i>Tritrichomonas foetus</i> while inhibiting fungi and bacteria
SPECIMEN	Rectal swab
TIME TO RESULT	Presumptive Positive: 15 - 30 minutes 2 - 4 days to positive
AVAILABILITY	5 pack 12-071-001 10 pack 12-071-002 100 pack 12-071-004



Learn more at
biomeddiagnostics.com

or call 866.964.6466
to speak to a representative