

# InPouch<sup>®</sup> TV

A highly selective culture media device for the detection of *Trichomonas vaginalis*, InPouch TV is a cost-effective, easy to use, broth media device that combines sample inoculation, incubation, and detection. Widely used in STI clinics, out-patient clinics, hospitals and remote, field-based healthcare for quick presumptive positives.

**PRINCIPLES**

- Direct inoculation to media at the POC
- Easy detection and identification

**COST SAVINGS**

- Conveniently collect samples, culture, transport and result in one device
- One-year shelf life at room temperature (18 - 25°C) from date of manufacture
- Oxygen resistant plastic pouch seals for culture and incubation

**SPECIFICITY**

- Specially formulated medium is selective for the growth of *T. vaginalis*, while inhibiting fungi and bacteria

**SPECIMEN**

Urine | Vaginal swab | CST | Seminal fluid

**RESULTS**

- Incubate at 35 - 37°C
- 15 - 30 minutes for presumptive positive

**CPT CODES**

87081 87088 87210

**AVAILABILITY**

10 pk		11-031-001
100 pk		11-031-002
Lab Kit (US)		11-035-001
Lab Kit (International)		11-035-002
TV Live Culture		11-041-003
TV Live Culture (10 pk)		11-051-004



Learn more at  
[biomeddiagnostics.com](http://biomeddiagnostics.com)

or call 866.964.6466  
 to speak to a representative