

Culture
Transport
Incubate
and Result
**IN ONE
DEVICE**

- InPouch can be viewed directly on the microscope stage
- Rapid turn-around time
- No subculturing required for detection
- FDA Clearance and CE IVD Mark on select products
- 12-month shelf life
- Reduces contamination risk and lab materials

Learn more at
biomeddiagnostics.com

Culture media
... **simplified**

Resealable Tabs

Prevents media dehydration
Small lab footprint

Durable Pouch

Oxygen-resistant
Designed for field use

**2 Chamber
Media Well**

Inoculate upper
chamber

Incubate & result in
lower

