



Certificate of Analysis
InTray™ COLOREX™ ESBL
Product Number: 11-173-001 and 11-173-002

Lot Number: 6SA054

Expiration Date: 05MAY2027

Test	Specification	Results
Appearance	Transparent, amber hue	Transparent, amber hue
Microbial load	No growth	No growth
pH	7.0 ± 0.2	7.0
Microbiology (ATCC®):	Colony Appearance	
<i>E. coli</i> BAA-196	Red	Red
<i>K. pneumoniae</i> 700603	Metallic blue	Metallic blue
<i>E. faecalis</i> 29212	Inhibited	Inhibited
<i>P. aeruginosa</i> 10145	Inhibited	Inhibited
<i>E. coli</i> 25922	Inhibited	Inhibited
<i>C. albicans</i> 60193	Inhibited	Inhibited
<i>S. aureus</i> 25923	Inhibited	Inhibited

QC Acceptance Date: 15MAY2026

Signature: *Megan Slama*