



**Certificate of Analysis**  
**InTray™ COLOREX™ GBS**

Product Number: 11-563-001, 11-563-002

Lot Number: 6SA055

Expiration Date: 06MAY2027

Test	Specification	Results
Appearance	Opalescent, clear	Opalescent, clear
Microbial load	No growth	No growth
pH	7.3 ± 0.2	7.22
<b>Microbiology (ATCC®):</b>	<b>Colony Appearance</b>	
<i>S. agalactiae</i> 12386	Growth, Mauve	Growth, Mauve
<i>S. agalactiae</i> 13813	Growth, Mauve	Growth, Mauve
<i>E. faecalis</i> 29212	Growth, Blue	Growth, Blue
<i>E. Coli</i> 25922	Inhibited	Inhibited
<i>C. albicans</i> 60193	Inhibited	Inhibited

**QC Acceptance Date:** 15MAY2026

**Signature:** *Megan Slama*