



Certificate of Analysis

Chocolate Eugon Agar with 10% Horse Blood

Product Number: 12-133-001, 12-133-002, 12-134-001, 12-134-002

Lot Number: 6SA059

Expiration Date: 11NOV2026

Test	Specification	Results
Appearance	Chocolate brown	Chocolate brown
Microbial load	No growth	No growth
pH	7.0 ± 0.1	6.91
Microbiology (STRAIN):	Colony Appearance	
<i>T. equigenitalis</i> 35865	Growth/oxichrome positive	Growth/oxichrome positive
<i>E. coli</i> 25922	Growth/oxichrome negative	Growth/oxichrome negative
<i>S. xylosus</i> 29971	Growth/oxichrome negative	Growth/oxichrome negative
<i>C. albicans</i> 90028	Growth/oxichrome variable	Growth/oxichrome positive

QC Acceptance Date: 18MAY2026

Signature: *Megan Slama*